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
Sophisticated Life Science Research Instrumentation


info@TSE-Systems.com


Transfer of validated behavioral paradigms into automated setup

Spatial learning and memory

A Morris water maze

B IntelliCage, Place Preference

Rsk2 Knockout Affects Emotional Behavior in the IntelliCage

Matthias Fischer¹ · Victoria Cabello¹ · Sandy Popp² · Sven Krackow³ ·
Leif Hommers¹ · Jürgen Deckert¹ · Klaus-Peter Lesch^{2,4,5} ·
Angelika G. Schmitt-Böhmer¹

Behav Genet

Behavioral tests in the IntelliCage

About 10–16 mice per IntelliCage (approximately equal numbers of KO and WT mice) were subjected to behavioral tests. Prior to each conditioning test, subjects were given an adaptation period of 6 days with all doors open and free access to the drinking bottles. Standard lab rodent food was provided ad lib throughout experiments. In addition, 2 days were given for the animals to adapt to nosepoking for door opening, i.e., doors were closed and opened for 7 s per visit only after the subject poked at a door (exception: 5 s for motor impulsivity).

Several experiments (all place learning schemes, anxiety test) were conducted using two drinking sessions during dark phase (21–23 and 3–5 h for place learning, 23–3 h for anxiety test) in order to increase learning motivation. Conditioning schemes were applied during drinking sessions only, IntelliCage doors stayed shut outside those sessions. In these cases, animals were given 2 days of adaptation to the schedule, i.e. doors opened at nosepoke for 7 s during drinking sessions, but kept close otherwise, before conditioning schemes were applied.

Activity and exploratory behavior

During the 6 days of free adaptation, the following behavioral measures were extracted for characterisation of spontaneous activity of each subject: NPvisits=Visits with Nosepoke without Licks/h; Lvisits=visits with licks/h; Svisits=visits without licks and nosepokes/h; NPVdur=median duration of visits w/np w/out lick; Nosepokes=mean number of nosepokes during visits w/np w/out licks; NPduration=median duration of such np during a visit; Licks=median number of licks per visit; Lduration=median duration of licking during a visit; Lcontact=median bottle cap contact time during a visit; Nocturnal=relative night activity [(visits during night – visits during day)/(total no. of visits)]; Repetitive=log(sum of observed returns to same corner/sum of expected such switches); Regularity=sqrt(sum of squared non-diag transition matrix residuals/sum of non-diag transition matrix observed values). We used the beginning of this 6 days adaptation period after placing the mice into the IntelliCage to investigate exploratory activity of a new environment by measuring time till first visit/nosepoke/lick, number of visits till first nosepoke/lick or time and number of visits before visiting, nosepoking or drinking in all four corners. Means are taken for frequencies with few levels when medians could strongly mislead in small samples and mask differences in large ones. Medians for durations and licks cut out outliers, like corner sleepover or spurious hits of sensors during exiting.

Behav Genet

Place learning

After free adaptation, nosepoke adaptation, and drinking session adaptation, a series of place learning experiments was carried out consecutively, with several experimental groups. In all cases, visits with nosepokes represented a trial that could either be correct or not (25% success by chance).

Place preference

Nosepokes opened the door only in one of the four corners for each subject for 6 sessions (3 days). Rewarded corner was randomly assigned to each subject and balanced between genotypes. After that, rewarded corner was swapped to the opposite corner (reversal) for another 6 sessions.

Place75

Consecutive to Place learning reversal, rewarded corner was set to the next corner in clockwise direction and reward given with a probability of 75%, for another 6 sessions.

Place change

For 6 sessions the rewarded corner was swapped randomly (excluding the prior rewarded corner), at every session.

Chaining

For 8 sessions, the rewarded corner was switched after each visit with nosepoke to the next corner in clockwise direction. Hence, subjects had to circle around the IntelliCage to get water during the drinking sessions. After that, direction of rewarded corner switch was reversed to counterclockwise (reversal) for another 16 sessions.

Patrolling

Same as chaining except that the rewarded corner switched to the next corner after each *correct* visit with nosepoke, i.e. only after a visit with nosepoke that led to door opening. Hence, in contrast to chaining, the position of the rewarding corner was not switched until it was successfully entered. After 12 such sessions (6 days), direction of rewarded corner switch was changed from clockwise to counterclockwise (reversal) for 16 sessions.

Motor impulsivity

During training phase (first 3 days) the door initially poked at with the first nosepoke of a visit opened 0.5, 1.5 or 2.5 s

after the initial nosepoke. The length of the delay varied randomly. After the delay 3 green LEDs were switched on on the side of the first nosepoke and the door opened at this side for 5 s. During the training phase (first 3 days) a nosepoke before the LEDs were switched on (premature nosepoke) had no effect.

In the following test phase (5 days) premature second nosepokes (before LED on) terminated the trial, i.e., the animals would have to leave the corner and visit again to start a new trial. Hence, animals had to learn, that after an initial nosepoke they had to wait with the second one till LED was switched on, not knowing if it would take 0.5, 1.5 or 2.5 s.

Impulsivity level was therefore measured by the proportion of correct responses (i.e., visits without premature nosepokes).

Cognitive impulsivity

In each corner one bottle was filled with water and one with 2% sucrose. A nosepoke at one door prevented opening of the other door during a visit. During 4 days mice were allowed to develop a preference for sucrose.

After that, the delay discounting condition scheme allowed drinking after nosepoking at the water bottle side without delay, as before, but at the sucrose side mice had to wait for some period of time till a nosepoke could open the door. The delay increased from 0.5 s every day by 0.5 s, up to 8 s delay. An earlier breakdown of preference for sucrose side would be taken as evidence of higher cognitive impulsivity in a group.

Chaining, patrolling and impulsivity were determined with similar designs as applied in Kobayashi et al. (2013).

Anxiety

In each corner one bottle was filled with water and one with 2% sucrose. WT and KO mice were assigned randomly and balanced within strains to two groups, respectively. One group had access only to the water bottles, the other group only to the sucrose bottles.

After general adaptation to the IntelliCage with water in all bottles and 3 days adaptation to the described sucrose and water side situation, eight days of airpuff-application followed. From 23 h till 3 h a nosepoke in one or the four corners initiated an airpuff for 1 s and with about 1 bar ("anxiety" phase). The anxiety corner was assigned randomly and balanced within strains to each mouse. During the intervening time periods no airpuffs were delivered ("inter" phase).

As an indicator of state anxiety, we measured the time it took a subject, after reception of the first airpuff, to nosepoke again in the airpuff delivering corner.

Cognitive Flexibility

frontiers in
BEHAVIORAL NEUROSCIENCE

ORIGINAL RESEARCH ARTICLE
published: 29 April 2014
doi: 10.3389/fnbeh.2014.00140

A novel automated behavioral test battery assessing cognitive rigidity in two genetic mouse models of autism

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
¹ Department of Neurophysiology, Nencki Institute of Experimental Biology, Warsaw, Poland

² Division of Functional Neuroanatomy, Institute of Anatomy, University of Zurich, Zurich, Switzerland

³ Department of Physiology, School of Laboratory Medicine, Kwazulu-Natal University, Durban, South Africa

	Adaptation phases		Training phases	
	Simple adaptation	Nosepoke adaptation	Place preference learning	Place reversal learning
Days	1–5	6–7	8–10	11–13
Liquid in bottles	Tap water	Tap water	10% sucrose	10% sucrose
Learning chambers with liquid access	All 4	All 4	One per animal	One per animal (opposite to the one accessible during the previous phase)
IntelliCage overview				

Effort choice behavior

Characterization of an alcohol addiction-prone phenotype in mice

Kasia Radwanska & Leszek Kaczmarek
Laboratory of Molecular Neurobiology, Nencki Institute, Poland


Figure 1 Experimental timeline